

Hemolysis Profile of a Novel Balloon-Filled Basket Pulsed Field Ablation Catheter

Jed A. Overmann¹, DVM, MS, DACVP, Mayara Marques², MD, Chris Lafean¹, BS, Catherine Pipenhagen¹, BS, Boyce L. Moon¹, MS, Atul Verma², MD

¹Abbott, St. Paul, MN, United States, ²McGill University Health Centre, Montreal, Quebec, Canada

Background

Pulsed field ablation (PFA) can cause hemolysis depending on the field delivered into the bloodstream. While any PFA system can theoretically cause hemolysis, the degree to which this happens is system dependent.




Figure 1. VOLT™ PFA Catheter, Sensor Enabled™

Objective

To assess the degree of hemolysis caused by a balloon-filled basket PFA catheter in an in vivo pre-clinical model.

Methods

Six canines were treated with a novel 12.5F, deflectable, variable diameter, balloon-filled, 8-spline-electrode basket PFA catheter (VOLT™ PFA Catheter, Sensor Enabled™, Abbott) (see Figure 1). Canines were survived for either 5-7 (N=3) or 30-32 days (N=3). All four pulmonary veins (PV) were treated in each animal with PFA dosing at twice the recommended upper limit (with respect to the number of applications) to create a challenge condition for hemolysis. For 20/24 PVs, this meant 16 applications per vein at 1800V. Six applications per vein of 1400V were delivered to the remaining 4 PVs (all RSPVs) due to proximity to the phrenic nerve. Samples were collected for plasma free hemoglobin (FH), hemoglobin (HGB) and serum creatinine (CR) at baseline (T0), after all PFA treatments in the left atrium (T1), and at termination at either 5-7 days or 30-32 days (TERM). Statistical analysis was completed using SigmaPlot 15.

Results

No acute complications were observed during PFA application, and all six canines survived until termination. T1 mean plasma FH concentration was statistically less than the clinical threshold of 50 mg/dL defining major hemolysis (one sample t-test, $p < 0.001$, with upper end of 95% confidence interval = 23.4 mg/dL)[1,2]. Plasma FH increased from a baseline (T0) mean of 6.6 ± 5.8 mg/dL to 16.4 ± 5.6 mg/dL immediately after all PFA applications in the left atrium (T1). This returned to 12.1 ± 6.5 mg/dL at TERM. For HGB, values were 13.0 ± 1.7 g/dL at T0, 11.5 ± 0.5 g/dL at T1, and 13.4 ± 1.6 g/dL at TERM. For CR, values were 0.9 ± 0.1 mg/dL, 0.8 ± 0.1 mg/dL, and 0.9 ± 0.1 mg/dL at T0, T1, and TERM respectively. The above results along with additional serum values are depicted in Table 1.

Table 1.

	T0	T1	TERM	p-value ³
Plasma FH (mg/dL) ^{1,2}	6.6±5.8	16.4±5.6	12.1±6.5	0.128
HGB (g/dL) ¹	13.0±1.7	11.5±0.5	13.4±1.6	0.014
			§	
Creatinine (mg/dL) ¹	0.9±0.1	0.8±0.1	0.9±0.1	0.126
Albumin (g/dL) ¹	2.9±0.3	2.5±0.1	3.0±0.2	<0.001
		*	§	
LDH (U/L) ¹	52±30	156±47	54±14	<0.001
		*	§	
BILI (mg/dL) ¹	0.1±0	0.1±0.1	0.1±0	0.402
CK (U/L) ¹	203±47	1362±537	467±693	<0.001
		*	§	

FH= free hemoglobin, HGB= hemoglobin, LDH = lactate dehydrogenase; BILI = bilirubin; CK = creatine kinase

¹ One Way Repeated Measures Analysis of Variance

² Plasma FH results only available for N=5 animals


³ p-value < 0.05 indicates the differences in mean values among timepoints are greater than would be expected by chance. Multiple comparison procedure (Bonferroni t-test) used to isolate timepoints that differed.

§ Statistically different ($p < 0.05$) vs. T1, based on Bonferroni t-test for pairwise comparisons

* Statistically different ($p < 0.05$) vs. T0, based on Bonferroni t-test for pairwise comparisons

Conclusion

Small amounts of hemolysis were noted after very high numbers of PFA applications using the VOLT™ PFA system. The amounts seen are less than those reported by some other PFA systems (see Figure 2).



Disclosures

Study was funded by Abbott. Authors Overmann, Lafean, Pipenhagen, and Moon are current Abbott employees. Dr. Verma has disclosures of advisory and grants from Biosense Webster, and Medtronic and advisory for Abbott and MedLumics. Dr. Marques has no disclosures.

Acknowledgements

Study was conducted at NAMSA, Minneapolis, MN

References

- [1] Dufour N, Radjou A, Thuong M. Hemolysis and Plasma Free Hemoglobin During Extracorporeal Membrane Oxygenation Support: From Clinical Implications to Laboratory Details. ASAIO J. 2020 Mar;66(3):239-246. doi: 10.1097/MAT.0000000000000974. PMID: 30985331.
- [2] ELSO Guidelines for Cardiopulmonary Extracorporeal Life Support. Extracorporeal Life Support Organization, Version 1.4 August 2017. Ann Arbor, MI, USA www.elseo.org
- [3] Nies M, Koruth JS, Mlček M, Watanabe K, Tibenská VC, Královce Š, Tejkl L, Neuzil P, Reddy VY. Hemolysis After Pulsed Field Ablation: Impact of Lesion Number and Catheter-Tissue Contact. Circ Arrhythm Electrophysiol. 2024 Jun;17(6):e012765. doi: 10.1161/CIRCEP.124.012765. Epub 2024 Apr 23. PMID: 38651357.
- [4] Mattison L, Verma A, Tarakji KG, Sigg DC. Hemolysis After Pulsed Field Ablation: The Role of Dose and Contact in an Acute Porcine Model. Circ Arrhythm Electrophysiol. 2024 Dec;17(12):e013317. doi: 10.1161/CIRCEP.124.013317. Epub 2024 Nov 27. PMID: 39601134.
- [5] De Smet MAJ, François C, De Becker B, Tavernier R, le Polain de Waroux JB, Knecht S, Duytschaever M. Intravascular haemolysis and acute kidney injury following atrial fibrillation ablation: a report using two different systems for pulsed field ablation. Europace. 2024 Oct 3;26(10):euae251. doi: 10.1093/europace/euae251. PMID: 39351957; PMCID: PMC11493093.